myUnity: What's New!

Release 2022.16




Release Dates


CY 2022 myUnity Senior Living Release Schedule (all dates are targets and subject to change)				
Release Value	GA Date (self-hosted)	UAT Install Date (Plexus-hosted)	Senior Living Prod Install Date (Plexus-hosted)	
2022.16	08/08	08/10	08/31	



Netsmart Connect

- NetsmartCONNECT
 Sign in
- Get Connected, join a Community Group.







Orders Census

Scheduling

Clinical

Financial

Personnel

Setup

Training Portal

Senior Living

MUE-13893

Overview

MUE-13893

Vitals Integration with Health Recovery Solutions (HRS)

The HL7 vitals integration with Health Recovery Solutions will allow the following data to be sent from Health Recovery Solutions to myUnity thru HL7 ORU messages.

Vital Signs:

- Systolic and diastolic blood pressure
- Oxygen saturation
- Heart rate
- Blood sugar
- Temperature
- Weight

Clinical notes

- Category Tele-monitoring
- SubCategory Notes


Contact your Client Alignment Executive for more information.

Senior Living

MUE-13893

Functionality MUE-13893 Vitals Integration with Health Recovery Solutions (HRS)


Telemonitoring device must be a staff assignment in Scheduling Assignments.






ribbon.

```
Admission: Pnehr-40869 SN 3/1/2020 - 7/26/2022 V
                                                                         3mo 🗸
                                                                                     4/26/2022 - 7/26/2022
                                         151/81 2 181.00 98.70 102.00 6 40.00
 Primary Payer: Amy Med B
```







Senior Living


MUE-11838

Functionality MUE-11838 Allow Lot Number to be entered on Administered Immunizations

For administered immunizations, only the lot number and observed results would be editable.









Senior Living


MUE-13755

Setup MUE-13755 Simple LTC Integration

Integration with SimpleLTC for MDS Scrubbing

Setup>General>Organization Parameters

- Select the Organization
- Select MDS/Assessment Scheduler
- MDS External Analysis Software
 - Select SimpleLTC





Senior Living


MUE-13755


Functionality MUE-13755 Simple LTC Integration


Integration with SimpleLTC for MDS Scrubbing


Clinical>MDS Assessments>

Submit to SimpleLTC will now appear












Senior Living

Overview

Regulatory News

Regulatory News from the NTST Regulatory Analyst has been placed on the top of the Release Notes page on the training portal. Latest three entries will appear on top and an archive page will exist for those earlier than the latest three.



Release Dates

CY 2022 myUnity Senior Living Release Schedule (all dates are targets and subject to change)				
Release Value	GA Date (self-hosted)	UAT Install Date (Plexus-hosted)	Senior Living Prod Install Date (Plexus-hosted)	
2022.16	08/08	08/10	08/31	

